

Human prostacycline, PGI2 ELISA Kit | PGI2 elisa kit

Human prostacycline, PGI2 ELISA Kit

CATALOG #
AAA15235

REACTIVITY
Human


PRODUCT SYNONYMS

prostacycline; PGI2; Human prostacycline; PGI2 ELISA Kit; PGI2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human PGI2. No significant cross-reactivity or interference between human PGI2 and analogues was observed.
Samples	Serum, plasma
Assay Type	Quantitative Competitive
Detection Range	50 pg/ml -2000 pg/ml
Sensitivity	Typically less than 50 pg/ml.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<15%. Three samples of known concentration were tested twenty times on one plate to assess
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<15%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGI2 elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with goat-anti-rabbit antibody. Standards or samples are added to the appropriate microtiter plate wells with an antibody specific for PGI2 and Horseradish Peroxidase (HRP) conjugated PGI2. The competitive inhibition reaction is launched between with HRP labeled PGI2 and unlabeled PGI2 with the antibody. A substrate solution is added to the wells and the color develops in opposite to the amount of PGI2 in the sample. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/prostacycline-pgi2/15235>