

TERF1 Interacting Nuclear Factor 2 (TINF2) Recombinant Protein | TINF2 recombinant protein

Recombinant TERF1 Interacting Nuclear Factor 2 (TINF2)

CATALOG # AAA153712	HOST E coli	APPLICATIONS Western Blot, SDS-Page
-------------------------------	-----------------------	---

PRODUCT SYNONYMS

TERF1 Interacting Nuclear Factor 2 (TINF2); Recombinant TERF1 Interacting Nuclear Factor 2 (TINF2); TERF1 Interacting Nuclear Factor 2; TINF2 recombinant protein (prokaryotic); Human; TIN2; TRF1-interacting nuclear protein 2; TINF2 recombinant protein

PRODUCT SPECIFICATIONS

Host	E coli
Purity/Purification	Greater than 95% by SDS-PAGE
Form/Format	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose. Freeze-dried powder
Sequence Positions	Ala2~Leu451
Applicable Applications for TINF2 recombinant protein	WB (Western Blot), SDS-PAGE, Immunogen, Positive Control
Expression System	Prokaryotic expression
Species	Human
Tag	N-terminal His
Endotoxin Level	<1.0EU per 1ug (determined by the LAL method)
Isoelectric Point	8
Subcellular Location	Nucleus, Chromosome
Research Area	Signal transduction
Preparation and Storage	Store at 2-8 degree C for one month, -80 degree C for 12 months Stable for 12 months if correctly stored.



ADDITIONAL INFORMATION

Product Categories/Family for TINF2 recombinant protein	Prokaryotic proteins
---	----------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	151101266
NCBI GeneID	26277
NCBI Accession #	NP_001092744.1
NCBI GenBank Nucleotide #	NM_001099274.1
UniProt Accession #	Q9BSI4; Q9H904; Q9UHC2
Molecular Weight	50,023 Da
NCBI Official Full Name	TERF1-interacting nuclear factor 2 isoform 1
NCBI Official Synonym Full Names	TERF1 (TRF1)-interacting nuclear factor 2
NCBI Official Symbol	TINF2
NCBI Official Synonym Symbols	TIN2; DKCA3
NCBI Protein Information	TERF1-interacting nuclear factor 2; TRF1-interacting nuclear protein 2
UniProt Protein Name	TERF1-interacting nuclear factor 2
UniProt Gene Name	TINF2
UniProt Synonym Gene Names	TIN2
UniProt Entry Name	TINF2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.