

Human alpha II-spectrin breakdown product SBDP145 ELISA Kit

Human alpha II-spectrin breakdown product SBDP145 ELISA Kit

CATALOG #
AAA15899

REACTIVITY
Human


PRODUCT SYNONYMS

alpha II-spectrin breakdown product SBDP145; Human alpha II-spectrin breakdown product SBDP145 ELISA Kit; Human alpha II-spectrin breakdown product SBDP145 ELISA kit; alpha II-spectrin breakdown product SBDP145 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human SBDP145. No significant cross-reactivity or interference between human SBDP145 and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cerebrospinal fluid (CSF), tissue homogenates
Detection Range	0.625 ng/ml-40 ng/ml
Sensitivity	< 0.156 ng/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for alpha II-spectrin breakdown product SBDP145 elisa kit

Principle of the Assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for SBDP145 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SBDP145 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SBDP145 is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SBDP145 bound in the initial step. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/alpha-ii-spectrin-breakdown-product-sbdp145/15899>