

Human Solube receptor for advanced glycation end products ELISA Kit | SRAGE elisa kit

Human Solube receptor for advanced glycation end products ELISA Kit

CATALOG #
AAA16123

REACTIVITY
Human


PRODUCT SYNONYMS

Solube receptor for advanced glycation end products; Human Solube receptor for advanced glycation end products ELISA Kit; SRAGE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of sRAGE. No significant cross-reactivity or interference between sRAGE and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between sRAGE and all the analogues, therefore, cross reaction may still exist in some cases.
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 pg/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SRAGE elisa kit	Intended Uses: This sRAGE ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human sRAGE. This ELISA kit for research use only, not for therapeutic or test applications! Principle of the Assay: sRAGE ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-sRAGE antibody and an sRAGE-HRP conjugate. The assay sample and buffer are incubated together with sRAGE-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the sRAGE concentration since sRAGE from samples and sRAGE-HRP conjugate compete for the anti-sRAGE antibody binding site. Since the number of sites is limited, as more sites are occupied by sRAGE from the sample, fewer sites are left to bind sRAGE-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The sRAGE concentration in each sample is interpolated from this standard curve.
Product Categories/Family for SRAGE elisa kit	Neurobiology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/srage-human-elisa-kits/solube-receptor-for-advanced-glycation-end-products/16123>