

Human deltaSIP-alpha ELISA Kit | deltaSIP-alpha elisa kit

Human deltaSIP-alpha (Delta Sleep Inducing Peptide alpha) ELISA Kit

CATALOG #
AAA164415

REACTIVITY
Human


PRODUCT SYNONYMS

deltaSIP-alpha; Human deltaSIP-alpha (Delta Sleep Inducing Peptide alpha) ELISA Kit; deltaSIP-alpha elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human deltaSIP-alpha. No significant cross-reactivity or interference between Human deltaSIP-alpha and analogues was observed. Note: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between Human deltaSIP-alpha and all the analogues, therefore, cross reaction may still exist.
Samples	Serum, Plasma, Biological Fluids
Assay Type	Competitive
Detection Range	125-8000pg/mL
Sensitivity	Min: 75pg/mL; Max: 8000pg/mL
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for deltaSIP-alpha elisa kit

Description: The kit is a competitive enzyme immunoassay for in vitro quantitative measurement of deltaSIP-alpha in human serum, plasma and other biological fluids

Principle of the Assay: This ELISA kit uses Competitive-ELISA as the method. The microtiter plate provided in this kit has been pre-coated with Human deltaSIP-alpha. During the reaction, Human deltaSIP-alpha in the sample or standard competes with a fixed amount of Human deltaSIP-alpha on the solid phase supporter for sites on the Biotinylated Detection Ab specific to Human deltaSIP-alpha. Excess conjugate and unbound sample or standard are washed from the plate, and Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then a TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The concentration of Human deltaSIP-alpha in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/deltasip-alpha/164415>