

Mouse NF-kappaB p65 ELISA Kit | NF-kappaB p65 elisa kit

Mouse NF-kappaB p65 (Nuclear Factor-kappaB p65) ELISA Kit

CATALOG #
AAA164612

REACTIVITY
Mouse


PRODUCT SYNONYMS

NF-kappaB p65; Mouse NF-kappaB p65 (Nuclear Factor-kappaB p65) ELISA Kit; NF-kappaB p65 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This kit recognizes Mouse NFE2L2 in samples. No significant cross-reactivity or interference between Mouse NFE2L2 and analogues was observed.
Sequence Length	969
Samples	Serum, Plasma And Other Biological Fluids
Assay Type	Quantitative Sandwich
Detection Range	0.31-20ng/mL
Sensitivity	0.19ng/mL.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Mouse NFE2L2 were tested 20 times on one plate, respectively
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Mouse NFE2L2 were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NF-kappaB p65 elisa kit	Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Mouse NFE2L2 concentrations in serum, plasma and other biological fluids. Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse NFE2L2. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse NFE2L2 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse NFE2L2, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of Mouse NFE2L2. You can calculate the concentration of Mouse NFE2L2 in the samples by comparing the OD of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	34577122
NCBI GeneID	4790
NCBI Accession #	NP_003989.2
NCBI GenBank Nucleotide #	NM_003998.3
UniProt Accession #	P19838; Q68D84; Q86V43; Q8N4X7; Q9NZC0; A8K5Y5; B3KVE8
Molecular Weight	85,520 Da
NCBI Official Full Name	nuclear factor NF-kappa-B p105 subunit isoform 1
NCBI Official Synonym Full Names	nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
NCBI Official Symbol	NFKB1
NCBI Official Synonym Symbols	p50; KBF1; p105; EBP-1; NF-kB1; NFKB-p50; NFKappaB; NF-kappaB; NFKB-p105; NF-kappa-B
NCBI Protein Information	nuclear factor NF-kappa-B p105 subunit; DNA binding factor KBF1; DNA-binding factor KBF1; NF-kappabeta; nuclear factor NF-kappa-B p50 subunit; nuclear factor kappa-B DNA binding subunit
UniProt Protein Name	Nuclear factor NF-kappa-B p105 subunit
UniProt Gene Name	NFKB1
UniProt Entry Name	NFKB1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.