

Human I-TAC ELISA Kit | I-TAC elisa kit

Human I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit

CATALOG #
AAA164967

REACTIVITY
Human


PRODUCT SYNONYMS

I-TAC; Human I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit; I-TAC elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human I-TAC. No significant cross-reactivity or interference between Human I-TAC and analogues was observed. Note: Limited by current skills and knowledge, it is impossible for us to complete the cross- reactivity detection between Human I-TAC and all the analogues, therefore, cross reaction may still exist.
Sequence Length	73
Samples	Serum, Plasma, Biological Fluids
Assay Type	Sandwich
Detection Range	62.5-4000pg/mL
Sensitivity	Min: 37.5pg/mL; Max: 4000pg/mL
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for I-TAC elisa kit	Description: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of I-TAC in human serum, plasma and other biological fluids Principle of the Assay: This ELISA kit uses Sandwich-ELISA as the method. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human I-TAC. Standards or samples are added to the appropriate micro ELISA plate wells and bound by the specific antibody. Then a biotinylated detection antibody specific for Human I-TAC and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human I-TAC, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of Human I-TAC. You can calculate the concentration of Human I-TAC in the samples by comparing the OD of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	14030443
NCBI GeneID	6373
UniProt Accession #	O14625; Q53YA3; Q92840
Molecular Weight	10,365 Da
NCBI Official Full Name	I-TAC, partial
NCBI Official Synonym Full Names	chemokine (C-X-C motif) ligand 11
NCBI Official Symbol	CXCL11
NCBI Official Synonym Symbols	IP9; H174; IP-9; b-R1; I-TAC; SCYB11; SCYB9B
NCBI Protein Information	C-X-C motif chemokine 11; beta-R1; interferon gamma-inducible protein 9; interferon-inducible T-cell alpha chemoattractant; small inducible cytokine B11; small inducible cytokine subfamily B (Cys-X-Cys), member 11; small inducible cytokine subfamily B (Cys-X-Cys), member 9B; small-inducible cytokine B11
UniProt Protein Name	C-X-C motif chemokine 11
UniProt Gene Name	CXCL11
UniProt Synonym Gene Names	ITAC; SCYB11; SCYB9B; IP-9; I-TAC
UniProt Entry Name	CXL11_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/i-tac/164967>