

Human TE ELISA Kit | TE elisa kit

Human TE (Telomerase) ELISA Kit

CATALOG #
AAA164975

REACTIVITY
Human


PRODUCT SYNONYMS

TE; Human TE (Telomerase) ELISA Kit; TE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes natural and recombinant Human TE. No significant cross-reactivity or interference between Human TE and analogues was observed.
Sequence Length	13
Samples	Serum, Plasma, Biological Fluids
Assay Type	Sandwich
Detection Range	0.156-10ng/mL
Sensitivity	Min: 0.094ng/mL; Max: 10ng/mL
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TE elisa kit	Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Human TE concentrations in serum, plasma and other biological fluids. Principle of the Assay This ELISA kit uses Sandwich-ELISA as the method. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to TE. Standards or samples are added to the appropriate micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for TE and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain TE, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of TE. You can calculate the concentration of TE in the samples by comparing the OD of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	3122066
NCBI Official Full Name	Testis ecdysiotropin peptide E
UniProt Protein Name	Testis ecdysiotropin peptide E
UniProt Gene Name	TE
UniProt Synonym Gene Names	TE
UniProt Entry Name	ECDE_LYMDI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/te/164975>