

Mouse Interleukin 35 ELISA Kit | IL-35 elisa kit

Mouse Interleukin 35 ELISA Kit

CATALOG # AAA16531	REACTIVITY Mouse
------------------------------	----------------------------


PRODUCT SYNONYMS

Interleukin 35; Mouse Interleukin 35 ELISA Kit; IL-35 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	50-1000pg/mL
Sensitivity	1.0pg/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IL-35 elisa kit	Intended Uses: This IL-35 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse IL-35. Principle of the Assay: IL-35 ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-IL-35 antibody and an IL-35-HRP conjugate. The assay sample and buffer are incubated together with IL-35-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the IL-35 concentration since IL-35 from samples and IL-35-HRP conjugate compete for the anti-IL-35 antibody binding site. Since the number of sites is limited, as more sites are occupied by IL-35 from the sample, fewer sites are left to bind IL-35-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The IL-35 concentration in each sample is interpolated from this standard curve.
Product Categories/Family for IL-35 elisa kit	Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/il-35-mouse-elisa-kits/interleukin-35/16531>