

Mouse VE-Cadherin ELISA Kit | VE-Cadherin elisa kit

Mouse VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit

CATALOG #
AAA165639

REACTIVITY
Mouse


PRODUCT SYNONYMS

VE-Cadherin; Mouse VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit; VE-Cadherin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This kit recognizes Mouse VE-Cadherin in samples. No significant cross-reactivity or interference between Mouse VE-Cadherin and analogues was observed.
Sequence Length	784
Samples	Serum, plasma and other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	0.16-10ng/mL
Sensitivity	0.1ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Mouse VE-Cadherin were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Mouse VE-Cadherin were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VE-Cadherin elisa kit	Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Mouse VE-Cadherin concentrations in serum, plasma and other biological fluids. Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse VE-Cadherin. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse VE-Cadherin and Avidin- Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse VE-Cadherin, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of Mouse VE-Cadherin. You can calculate the concentration of Mouse VE-Cadherin in the samples by comparing the OD of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	599834
NCBI GeneID	1003
UniProt Accession #	P33151; Q4VAI5; Q4VAI6
Molecular Weight	74,499 Da
NCBI Official Full Name	VE-cadherin
NCBI Official Synonym Full Names	cadherin 5, type 2 (vascular endothelium)
NCBI Official Symbol	CDH5
NCBI Official Synonym Symbols	7B4; CD144
NCBI Protein Information	cadherin-5; 7B4 antigen; VE-cadherin; cadherin 5, type 2, VE-cadherin (vascular epithelium); cd144 antigen; endothelial-specific cadherin; vascular endothelial cadherin
UniProt Protein Name	Cadherin-5
UniProt Gene Name	CDH5
UniProt Synonym Gene Names	VE-cadherin
UniProt Entry Name	CADH5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.