

Human Phosphorylation P38 Mitogen activated Protein Kinase ELISA Kit | P P38MAPK elisa kit

Human Phosphorylation P38 Mitogen activated Protein Kinase ELISA Kit

CATALOG # AAA16598	REACTIVITY Human
------------------------------	----------------------------


PRODUCT SYNONYMS

Phosphorylation P38 Mitogen activated Protein Kinase; Human Phosphorylation P38 Mitogen activated Protein Kinase ELISA Kit; P P38MAPK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive or Sandwich
Detection Range	50-1000ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for P P38MAPK elisa kit	Intended Uses: This PMAPKP38 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human PMAPKP38. This ELISA kit for research use only! Principle of the Assay: PMAPKP38 ELISA kit applies the quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with a monoclonal antibody specific for PMAPKP38. Standards or samples are then added to the microtiter plate wells and PMAPKP38 if present, will bind to the antibody pre-coated wells. In order to quantitatively determine the amount of PMAPKP38 present in the sample, a standardized preparation of horseradish peroxidase (HRP)-conjugated polyclonal antibody, specific for PMAPKP38 are added to each well to "sandwich" the PMAPKP38 immobilized on the plate. The microtiter plate undergoes incubation, and then the wells are thoroughly washed to remove all unbound components. Next, substrate solutions are added to each well. The enzyme (HRP) and substrate are allowed to react over a short incubation period. Only those wells that contain PMAPKP38 and enzyme-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The PMAPKP38 concentration in each sample is interpolated from this standard curve.
Product Categories/Family for P P38MAPK elisa kit	Cell Biology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/p-p38mapk-human-elisa-kits/phosphorylation-p38-mitogen-activated-protein-kinase/16598>