

Human Fcgamma (Fc Fragment of IgG) CLIA Kit | Fcg clia kit

Human Fcgamma (Fc Fragment of IgG) CLIA Kit

CATALOG #
AAA170875

REACTIVITY
Human


PRODUCT SYNONYMS

Fcgamma (Fc Fragment of IgG); Human Fcgamma (Fc Fragment of IgG) CLIA Kit; Human Fcg (Fc Fragment of IgG) CLIA Kit; Human Fcg; Fcg clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes natural and recombinant Fcgamma. No significant cross-reactivity or interference between Fcgamma and analogues was observed.
Assay Type	Sandwich
Detection Range	62.5-4000pg/mL
Sensitivity	Min: 37.5pg/mL; Max: 4000pg/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Fcg clia kit

Intended Uses: This CLIA kit applies to the in vitro quantitative determination of Fcgamma concentrations in serum, plasma and other biological fluids.

Principle of the Assay: This kit uses Sandwich-CLIA as the method. The micro CLIA plate provided in this kit has been pre-coated with an antibody specific to Fcgamma. Standards or samples are added to the appropriate micro CLIA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Fcgamma and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Fcgamma, biotinylated detection antibody and Avidin-HRP conjugate will appear fluorescence. The Relative light unit (RLU) value is measured spectrophotometrically by the Chemiluminescence immunoassay analyzer. The RLU value is positively associated with the concentration of Fcgamma. You can calculate the concentration of Fcgamma in the samples by comparing the RLU value of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

