

Human PGR (Progesterone Receptor) CLIA Kit | PGR clia kit

Human PGR (Progesterone Receptor) CLIA Kit

CATALOG #
AAA171091

REACTIVITY
Human


PRODUCT SYNONYMS

PGR (Progesterone Receptor); Human PGR (Progesterone Receptor) CLIA Kit; PGR clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes natural and recombinant PGR. No significant cross-reactivity or interference between PGR and analogues was observed.
Sequence Length	913
Assay Type	Sandwich
Detection Range	62.5-4000pg/mL
Sensitivity	Min: 37.5pg/mL; Max: 4000pg/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGR clia kit

Intended Uses: This CLIA kit applies to the in vitro quantitative determination of PGR concentrations in serum, plasma and other biological fluids.

Principle of the Assay: This kit uses Sandwich-CLIA as the method. The micro CLIA plate provided in this kit has been pre-coated with an antibody specific to PGR. Standards or samples are added to the appropriate micro CLIA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for PGR and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain PGR, biotinylated detection antibody and Avidin-HRP conjugate will appear fluorescence. The Relative light unit (RLU) value is measured spectrophotometrically by the Chemiluminescence immunoassay analyzer. The RLU value is positively associated with the concentration of PGR. You can calculate the concentration of PGR in the samples by comparing the RLU value of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	329112574
NCBI GeneID	100848466
NCBI Accession #	NP_001192285.1
NCBI GenBank Nucleotide #	NM_001205356.1
NCBI Official Full Name	progesterone receptor
NCBI Official Symbol	PGR
NCBI Protein Information	progesterone receptor

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/pgr-progesterone-receptor/171091>