

Human COL1 (Collagen Type I) CLIA Kit | COL1 clia kit

Human COL1 (Collagen Type I) CLIA Kit

CATALOG #
AAA171379

REACTIVITY
Human


PRODUCT SYNONYMS

COL1 (Collagen Type I); Human COL1 (Collagen Type I) CLIA Kit; COL1; COL1 clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	885
Detection Range	0.625~40ng/mL
Sensitivity	0.375ng/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for COL1 clia kit

The micro CLIA plate provided in this kit has been pre-coated with an antibody specific to Human COL1. Standards or samples are added to the appropriate micro CLIA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human COL1 and Avidin-Horseradish Peroxidase (HRP) conjugate is added to each micro plate well successively and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human COL1, biotinylated detection antibody and Avidin-HRP conjugate will appear fluorescence. The Relative light unit(RLU) value is measured spectrophotometrically by the Chemiluminescence immunoassay analyzer. The RLU value is positively associated with the concentration of Human COL1. You can calculate the concentration of Human COL1 in the samples by comparing the RLU value of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	119615037
NCBI GeneID	1277
Molecular Weight	138,941 Da
NCBI Official Full Name	collagen, type I, alpha 1, isoform CRA_a

NCBI Official Synonym Full Names	collagen, type I, alpha 1
NCBI Official Symbol	COL1A1
NCBI Official Synonym Symbols	OI1; OI2; OI3; OI4; EDSC
NCBI Protein Information	collagen alpha-1(I) chain
UniProt Protein Name	Collagen alpha-1(I) chain
UniProt Gene Name	COL1A1
UniProt Entry Name	CO1A1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/col1-collagen-type-i/171379>