

## Human Vasodilator Stimulated Phosphoprotein ELISA Kit | VASP elisa kit

Human VASP (Vasodilator Stimulated Phosphoprotein) ELISA kit

CATALOG #  
**AAA171414**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Vasodilator Stimulated Phosphoprotein; Human VASP (Vasodilator Stimulated Phosphoprotein) ELISA kit; VASP elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                                           |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                                                                                     |
| Specificity             | This kit recognizes Human VASP in samples. No significant cross-reactivity or interference between Human VASP and analogues was observed.                                 |
| Sequence Length         | 380                                                                                                                                                                       |
| Assay Type              | Quantitative Sandwich                                                                                                                                                     |
| Samples                 | Serum, plasma and other biological fluids                                                                                                                                 |
| Detection Range         | 78.13-5000pg/mL                                                                                                                                                           |
| Sensitivity             | 46.88pg/mL                                                                                                                                                                |
| Intra-assay Precision   | Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Human VASP were tested 20 times on one plate, respectively.               |
| Inter-assay Precision   | Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Human VASP were tested on 3 different plates, 20 replicates in each plate. |
| Preparation and Storage | Store at 4 degree C.                                                                                                                                                      |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for VASP elisa kit | <p>Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human VASP. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human VASP and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human VASP, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Human VASP. You can calculate the concentration of Human VASP in the samples by comparing the OD of the samples to the standard curve.</p> <p>Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Human VASP concentrations in serum, plasma and other biological fluids.</p> |
|------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

NCBI AND UNIPROT INFORMATION

|                                  |                                       |
|----------------------------------|---------------------------------------|
| NCBI GI #                        | 4507869                               |
| NCBI GeneID                      | 7408                                  |
| NCBI Accession #                 | NP_003361.1                           |
| NCBI GenBank Nucleotide #        | NM_003370.3                           |
| Molecular Weight                 | 39,830 Da                             |
| NCBI Official Full Name          | vasodilator-stimulated phosphoprotein |
| NCBI Official Synonym Full Names | vasodilator-stimulated phosphoprotein |
| NCBI Official Symbol             | VASP                                  |
| NCBI Protein Information         | vasodilator-stimulated phosphoprotein |
| UniProt Protein Name             | Vasodilator-stimulated phosphoprotein |
| UniProt Gene Name                | VASP                                  |
| UniProt Synonym Gene Names       | VASP                                  |
| UniProt Entry Name               | VASP_HUMAN                            |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/vasp-human-elisa-kits/vasodilator-stimulated-phosphoprotein/171414>