

Super Plus High Sensitive and Rapid Immunohistochemical Kit (pH6.0) IHC Kit

Super Plus High Sensitive and Rapid Immunohistochemical Kit (pH6.0)

CATALOG #
AAA177766

PRODUCT SYNONYMS

Super Plus High Sensitive and Rapid Immunohistochemical Kit (pH6.0); Super Plus High Sensitive and Rapid Immunohistochemical Kit (pH6.0) ihc kit

PRODUCT SPECIFICATIONS

Preparation and Storage

Store at 2~8 degree C, away from light. Avoid of freezing. Valid for 12 months. The reagents are valid within 6 months after opening.

ADDITIONAL INFORMATION

Related Product Information for Super Plus High Sensitive and Rapid Immunohistochemical Kit (pH6.0) ihc kit

Background: Dewaxing/Antigen Retrieval Buffer (pH6.0) (20×) [SP Reagent 6A] is a 20× concentrated solution. Dilute with DI water to 1 × Dewaxing/Antigen Retrieval working solution before use. For example: Take 10 mL Dewaxing/Antigen Retrieval Buffer (pH6.0) (20×), dilute with DI water to 200 mL.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/ihc-kits/super-plus-high-sensitive-and-rapid-immunohistochemical-kit-ph6-0/177766>