

Human IgA (Immunoglobulin A) ELISA Kit | IgA elisa kit

QuickKey Pro Human IgA (Immunoglobulin A) ELISA Kit

CATALOG #
AAA177884

REACTIVITY
Human


PRODUCT SYNONYMS

IgA (Immunoglobulin A); QuickKey Pro Human IgA (Immunoglobulin A) ELISA Kit; Immunoglobulin A; IgA; IgA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes Human IgA in samples. No significant cross-reactivity or interference between Human IgA and analogues was observed
Assay Type	Quantitative Sandwich
Samples	Human serum, plasma, saliva and urine samples
Detection Range	0.78-50 ng/mL
Sensitivity	0.14 ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Human IgA were tested 20 times on one plate, respectively.
Preparation and Storage	An unopened kit can be stored at 2-8° for 6 months. After test, the unused wells and reagents should be stored according to the table/manual.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IgA elisa kit	Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human IgA. Samples (or Standards) and Horseradish Peroxidase (HRP) linked antibody specific for Human IgA are added to the micro ELISA plate wells. Human IgA in samples (or standards) combines with the coated antibody and HRP linked detection antibody special to Human IgA. Excess conjugate and unbound sample or standard are washed from the plate. The substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450±2 nm. The OD value is proportional to the concentration of Human IgA. The concentration of Human IgA in the samples is then determined by comparing the OD of the samples to the standard curve. Intended Uses: This ELISA kit applies to the in vitro quantitative determination of IgA concentrations in Human serum, plasma, saliva and urine samples. Please contact tech-support for other sample type detection..
---	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/iga-immunoglobulin-a/177884>