

Rabbit Pg (Progesterone) ELISA Kit | Pg elisa kit

QuickKey Pro Rabbit Pg (Progesterone) ELISA Kit

CATALOG #
AAA177907

REACTIVITY
Rabbit


PRODUCT SYNONYMS

Pg (Progesterone); QuickKey Pro Rabbit Pg (Progesterone) ELISA Kit; Progesterone; Pg; Pg elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Specificity	This kit recognizes Pg in Rabbit samples. No significant cross-reactivity or interference between Rabbit Pg and analogues was observed
Assay Type	Quantitative Competitive
Samples	Rabbit serum and plasma samples
Detection Range	0.16-10 ng/mL
Sensitivity	0.05 ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Rabbit Pg were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Rabbit Pg were tested on 3 different plates, 20 replicates in each plate, respectively.
Preparation and Storage	An unopened kit can be stored at 2-8° for 6 months. After test, the unused wells and reagents should be stored according to the table/manual.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Pg elisa kit	Principle of the Assay: This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Rabbit Pg. Samples (or Standards) and Horseradish Peroxidase (HRP) linked antibody specific for Rabbit Pg are added to the micro ELISA plate wells. Rabbit Pg in samples (or standards) competes with a fixed amount of Pg on the solid phase supporter for sites on the HRP linked detection antibody specific to Pg. Excess conjugate and unbound sample or standard are washed from the plate. The substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450±2 nm. The concentration of Rabbit Pg in the samples is then determined by comparing the OD of the samples to the standard curve. Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Pg concentrations in Rabbit serum and plasma samples. Please contact tech-support for other sample type detection.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/pg-progesterone/177907>