

Fluor 680 Labeling Kit

Fluor 680 Labeling Kit

CATALOG #
AAA178189

PRODUCT SYNONYMS


Fluor 680 Labeling Kit; Fluor 680 Labeling Kit (30KD ~ 100KD protein labeling); Fluor 680 Labeling Kit labeling kit

PRODUCT SPECIFICATIONS

Preparation and Storage

2~8 degree C, shading light

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Fluor 680 Labeling Kit labeling kit

Background: The MBS Fluor 680 Labeling Kit provides all the reagents required for labeling, which can label proteins containing primary amino-group ((-NH₂) molecules simply and effectively.

Principle of the Assay: Within a certain pH range, MBS # Fluor 680 specifically reacts with primary amino groups (N-terminal and lysine residue side chains) to form a stable amide bond, so as to realize the coupling of MBS # Fluor 680 with protein.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/labeling-kits/fluor-680-labeling-kit/178189>