

Tartrate Resistant Acid Phosphatase (TRAP) Assay Kit | TRAP assay kit

Tartrate Resistant Acid Phosphatase (TRAP) Activity Assay Kit

CATALOG #
AAA178631

PRODUCT SYNONYMS

Tartrate Resistant Acid Phosphatase (TRAP); Tartrate Resistant Acid Phosphatase (TRAP) Activity Assay Kit; TRAP assay kit

PRODUCT SPECIFICATIONS

Samples	Serum, plasma, urine, animal tissue and cell samples
Detection Range	9.22- 100 U/L
Sensitivity	9.22 U/L
Intra-assay Precision	2
Inter-assay Precision	2.3
Preparation and Storage	-20 degree C

ADDITIONAL INFORMATION

Related Product Information for TRAP assay kit	Principle of the Assay: The chromogenic agent can be catalyzed by acid phosphatase to produce p-nitrophenol in acidic condition. P-nitrophenol has a maximum absorption at 405 nm. The activity of TRAP can be calculated by measur.iDg the amount of produced p-nitrophenol. The activity of acid phosphatase being detected in the presence of tartaric acid is regarded as TRAP activity. Intended Uses: This kit can measure tartrate resistant acid phosphatase (TRAP) activity in serum, plasma, urine, animal tissue and cell samples.
Product Categories/Family for TRAP assay kit	Metabolic Diseases

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/assay-kits/tartrate-resistant-acid-phosphatase-trap/178631>