

Rabbit Folic acid (FA) ELISA Kit | FA elisa kit

Rabbit Folic acid (FA) ELISA Kit

CATALOG # AAA179823	REACTIVITY Rabbit
-------------------------------	-----------------------------


PRODUCT SYNONYMS

Folic acid (FA); Rabbit Folic acid (FA) ELISA Kit; FA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Specificity	No cross-reaction with other factors.
Sequence Length	341
Samples	Serum, plasma or cell culture supernatant and organizations in the natural and recombinant FA concentration.
Assay Type	Sandwich
Detection Range	20 ng/ml-0.312 ng/ml
Sensitivity	0.06 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FA elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Rabbit FA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	334184341
NCBI GeneID	816615

NCBI Accession #	NP_001189564.1
NCBI GenBank Nucleotide #	NM_001202635.1
Molecular Weight	32,137 Da
NCBI Official Full Name	folic acid binding / transferase
NCBI Official Symbol	AT2G20830
NCBI Official Synonym Symbols	F5H14.20; F5H14_20
NCBI Protein Information	folic acid binding / transferase
UniProt Protein Name	Expressed protein
UniProt Entry Name	Q9SKT4_ARATH

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/folic-acid-fa/179823>