

Rabbit Decorin (DCN) ELISA Kit | DCN elisa kit

Rabbit Decorin (DCN) ELISA Kit

CATALOG #
AAA180084

REACTIVITY
Rabbit


PRODUCT SYNONYMS

Decorin (DCN); Rabbit Decorin (DCN) ELISA Kit; DCN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Sequence Length	354
Assay Type	Sandwich
Detection Range	100 ng/mL-1.56 ng/mL
Sensitivity	Up to 0.5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DCN elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Rabbit DCN monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	54020664
NCBI GeneID	29139
NCBI Accession #	NP_077043.1

NCBI GenBank Nucleotide #	NM_024129.1
Molecular Weight	39,805 Da
NCBI Official Full Name	decorin
NCBI Official Synonym Full Names	decorin
NCBI Official Symbol	Dcn
NCBI Protein Information	decorin; DSPG; PG40; PG-S2; bone proteoglycan II; dermatan sulfate proteoglycan-II (decorin)
UniProt Protein Name	Decorin
UniProt Gene Name	Dcn
UniProt Synonym Gene Names	DSPG
UniProt Entry Name	PGS2_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/decorin-dcn/180084>