

Rabbit Anti-Insulin Antibody (anti-Ins) ELISA Kit | Ins elisa kit

Rabbit Anti-Insulin Antibody (anti-Ins) ELISA Kit

CATALOG #
AAA180170

REACTIVITY
Rabbit


PRODUCT SYNONYMS

Anti-Insulin Antibody (anti-Ins); Rabbit Anti-Insulin Antibody (anti-Ins) ELISA Kit; Ins elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Specificity	No cross-reaction with other factors.
Sequence Length	123
Samples	Rabbit serum, plasma or cell culture supernatant and organizations in the natural and recombinant anti-Ins concentration.
Assay Type	Sandwich
Detection Range	200 U/ml-3.12 U/ml
Sensitivity	0.6 U/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Ins elisa kit

Intended Uses: This product is suitable for in vitro quantitative detection of Rabbit serum, plasma or cell culture supernatant and organizations in the natural and recombinant anti-Ins concentration.

Principle of the Assay||This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Rabbit anti-Ins monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	425688
NCBI Official Full Name	anti-insulin 193020 monoclonal antibody heavy chain variable region, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/anti-insulin-antibody-anti-ins/180170>