

Rabbit Monokine Induced By Interferon Gamma (Mlg) ELISA Kit | Mlg elisa kit

Rabbit Monokine Induced By Interferon Gamma (Mlg) ELISA Kit

CATALOG # AAA180193	REACTIVITY Rabbit
-------------------------------	-----------------------------


PRODUCT SYNONYMS

Monokine Induced By Interferon Gamma (Mlg); Rabbit Monokine Induced By Interferon Gamma (Mlg) ELISA Kit; Mlg elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Specificity	No cross-reaction with other factors.
Sequence Length	125
Samples	Rabbit serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	2000 pg/ml-31.2 pg/ml
Sensitivity	Up to 12 pg/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Mlg elisa kit	Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Rabbit Mlg monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	161579144
NCBI GeneID	513990

Molecular Weight	14,106 Da
NCBI Official Full Name	monokine induced by interferon gamma
NCBI Official Synonym Full Names	chemokine (C-X-C motif) ligand 9<
NCBI Official Symbol	CXCL9
NCBI Protein Information	C-X-C motif chemokine 9; small-inducible cytokine B9; monokine induced by interferon gamma
UniProt Protein Name	C-X-C motif chemokine 9
UniProt Gene Name	CXCL9
UniProt Entry Name	CXCL9_BOVIN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/monokine-induced-by-interferon-gamma-mig/180193>