

Human SET domain containing 5 ELISA Kit | SETD5 elisa kit

Human SET domain-containing protein 5, SETD5 ELISA Kit

CATALOG #
AAA180403

REACTIVITY
Human


PRODUCT SYNONYMS

SET domain containing 5; Human SET domain-containing protein 5; SETD5 ELISA Kit; SETD5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Samples	Human serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	10 ng/ml-0.156 ng/ml
Sensitivity	Up to 0.05 ng/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SETD5 elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human SETD5 monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	122937492
NCBI GeneID	55209
NCBI Accession #	NP_001073986.1
NCBI GenBank Nucleotide #	NM_001080517.2

UniProt Accession #	Q9C0A6; Q6AI17; Q8WUB6; Q9H3X4; Q9H6V7; Q9H7S3; Q9NVI9
Molecular Weight	157,515 Da
NCBI Official Full Name	SET domain-containing protein 5 isoform 1
NCBI Official Synonym Full Names	SET domain containing 5
NCBI Official Symbol	SETD5
NCBI Protein Information	SET domain-containing protein 5
UniProt Protein Name	SET domain-containing protein 5
UniProt Gene Name	SETD5
UniProt Synonym Gene Names	KIAA1757
UniProt Entry Name	SETD5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/set-domain-containing-5/180403>