

Cat butyrylcholinesterase ELISA Kit | BCHE elisa kit

Cat Cholinesterase, BCHE ELISA Kit

CATALOG #
AAA180456

REACTIVITY
Cat


PRODUCT SYNONYMS

butyrylcholinesterase; Cat Cholinesterase; BCHE ELISA Kit; BCHE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Cat
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BCHE elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Feline BCHE monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	57163907
NCBI GeneID	493960
NCBI Accession #	NP_001009364.1

NCBI GenBank Nucleotide #	NM_001009364.1
UniProt Accession #	O62760
Molecular Weight	68,328 Da
NCBI Official Full Name	cholinesterase
NCBI Official Symbol	BCHE
NCBI Protein Information	cholinesterase; choline esterase II; pseudocholinesterase; butyrylcholine esterase; acylcholine acylhydrolase
UniProt Protein Name	Cholinesterase
UniProt Gene Name	BCHE
UniProt Entry Name	CHLE_FELCA

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/cat-elisa-kits/butyrylcholinesterase/180456>