

Mouse Vasodilator-stimulated phosphoprotein ELISA Kit | VASP elisa kit

Mouse Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit

CATALOG #
AAA180524

REACTIVITY
Mouse


PRODUCT SYNONYMS

Vasodilator-stimulated phosphoprotein; Mouse Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit; VASP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	No cross-reaction with other factors.
Sequence Length	384
Samples	Mouse Serum, Plasma Or Cell Culture Supernatant And Organizations In The Natural And Recombinant Vaca Concentration
Assay Type	Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	0.5 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VASP elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Mouse VacA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	50979054
NCBI GeneID	403936

NCBI Accession #	NP_001003256.1
NCBI GenBank Nucleotide #	NM_001003256.1
Molecular Weight	40,413 Da
NCBI Official Full Name	vasodilator-stimulated phosphoprotein
NCBI Official Symbol	VASP
NCBI Protein Information	vasodilator-stimulated phosphoprotein
UniProt Protein Name	Vasodilator-stimulated phosphoprotein
UniProt Gene Name	VASP
UniProt Synonym Gene Names	VASP
UniProt Entry Name	VASP_CANLF

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/vasp-mouse-elisa-kits/vasodilator-stimulated-phosphoprotein/180524>