

Canine Gelsolin ELISA Kit | GS elisa kit

Canine Gelsolin (GS) ELISA Kit

CATALOG #
AAA180542

REACTIVITY
Canine


PRODUCT SYNONYMS

Gelsolin; Canine Gelsolin (GS) ELISA Kit; GS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Sequence Length	780
Samples	Canine Serum, Plasma Or Cell Culture Supernatant And Organizations In The Natural And Recombinant Grbeta Concentration
Assay Type	Sandwich
Detection Range	20 ng/ml-0.312 ng/ml
Sensitivity	0.06 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GS elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Canine GRbeta monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

51854227

NCBI GeneID	296654
NCBI Accession #	NP_001004080.1
NCBI GenBank Nucleotide #	NM_001004080.1
Molecular Weight	80,929 Da
NCBI Official Full Name	gelsolin
NCBI Official Synonym Full Names	gelsolin
NCBI Official Symbol	Gsn
NCBI Protein Information	gelsolin
UniProt Protein Name	Gelsolin
UniProt Gene Name	Gsn
UniProt Synonym Gene Names	ADF
UniProt Entry Name	GELS_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/gs-canine-elisa-kits/gelsolin/180542>