

Canine P-Selectin ELISA Kit | P-Sel elisa kit

Canine P-Selectin (P-Sel) ELISA Kit

CATALOG #
AAA180559

REACTIVITY
Canine


PRODUCT SYNONYMS

P-Selectin; Canine P-Selectin (P-Sel) ELISA Kit; P-Sel elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	830
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



ADDITIONAL INFORMATION

Related Product Information for P-Sel elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine P-Sel monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	157419154
NCBI GeneID	6403
NCBI Accession #	NP_002996.2

NCBI GenBank Nucleotide #	NM_003005.3
Molecular Weight	90,834 Da
NCBI Official Full Name	P-selectin
NCBI Official Synonym Full Names	selectin P
NCBI Official Symbol	SELP
NCBI Official Synonym Symbols	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
NCBI Protein Information	P-selectin
UniProt Protein Name	P-selectin
UniProt Gene Name	SELP
UniProt Synonym Gene Names	GMRP; GRMP; GMP-140; LECAM3; PADGEM
UniProt Entry Name	LYAM3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/p-sel-canine-elisa-kits/p-selectin/180559>