

Rat Pro Brain Natriuretic Peptide ELISA Kit | ProBNP elisa kit

Rat Pro Brain Natriuretic Peptide (ProBNP) ELISA Kit

CATALOG #
AAA180600

REACTIVITY
Rat


PRODUCT SYNONYMS

Pro Brain Natriuretic Peptide; Rat Pro Brain Natriuretic Peptide (ProBNP) ELISA Kit; ProBNP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Specificity	No cross-reaction with other factors.
Sequence Length	760
Assay Type	Sandwich
Detection Range	2000 pg/mL-31.2 pg/mL
Sensitivity	Up to 12 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ProBNP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Rat ProBNP monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	292495099
NCBI GeneID	2191

Molecular Weight	26,954 Da
NCBI Official Full Name	Prolyl endopeptidase FAP
NCBI Official Synonym Full Names	fibroblast activation protein alpha
NCBI Official Symbol	FAP
NCBI Official Synonym Symbols	FAPA; SIMP; DPPIV
NCBI Protein Information	prolyl endopeptidase FAP
UniProt Protein Name	Prolyl endopeptidase FAP
UniProt Gene Name	FAP
UniProt Synonym Gene Names	APCE
UniProt Entry Name	SEPR_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/probnp-rat-elisa-kits/pro-brain-natriuretic-peptide/180600>