

Canine Phospholipase C ELISA Kit | PLC elisa kit

Canine Phospholipase C (PLC) ELISA Kit

CATALOG #
AAA180608

REACTIVITY
Canine


PRODUCT SYNONYMS

Phospholipase C; Canine Phospholipase C (PLC) ELISA Kit; PLC elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	1252
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PLC elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine PLC monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	190036
NCBI GeneID	5336
Molecular Weight	147,870 Da

NCBI Official Full Name	phospholipase C
NCBI Official Synonym Full Names	phospholipase C gamma 2
NCBI Official Symbol	PLCG2
NCBI Official Synonym Symbols	FCAS3; APLAID; PLC-IV; PLC-gamma-2
NCBI Protein Information	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-2
UniProt Protein Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-2
UniProt Gene Name	PLCG2
UniProt Synonym Gene Names	PLC-IV; PLC-gamma-2
UniProt Entry Name	PLCG2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/plc-canine-elisa-kits/phospholipase-c/180608>