

Canine Prostacyclin ELISA Kit | PGI2 elisa kit

Canine Prostacyclin (PGI2) ELISA Kit

CATALOG #
AAA180614

REACTIVITY
Canine


PRODUCT SYNONYMS

Prostacyclin; Canine Prostacyclin (PGI2) ELISA Kit; PGI2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	386
Assay Type	Sandwich
Detection Range	2000 pg/mL-31.2 pg/mL
Sensitivity	Up to 12 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGI2 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine PGI2 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4506263
NCBI GeneID	5739

NCBI Accession #	NP_000951.1
NCBI GenBank Nucleotide #	NM_000960.3
Molecular Weight	40,956 Da
NCBI Official Full Name	prostacyclin receptor
NCBI Official Synonym Full Names	prostaglandin I2 (prostacyclin) receptor (IP)
NCBI Official Symbol	PTGIR
NCBI Official Synonym Symbols	IP; PRIPR
NCBI Protein Information	prostacyclin receptor
UniProt Protein Name	Prostacyclin receptor
UniProt Gene Name	PTGIR
UniProt Synonym Gene Names	PRIPR; PGI receptor; PGI2 receptor
UniProt Entry Name	PI2R_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/pgi2-canine-elisa-kits/prostacyclin/180614>