

Canine Amylin ELISA Kit | IAPP elisa kit

Canine Amylin ELISA Kit

CATALOG #
AAA180758

REACTIVITY
Canine


PRODUCT SYNONYMS

Amylin; Canine Amylin ELISA Kit; IAPP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	62
Samples	Canine serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	100 ng/mL-1.56 ng/mL
Sensitivity	Up to 0.5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IAPP elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine Amylin monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	457131
NCBI GeneID	3375

Molecular Weight	9,806 Da
NCBI Official Full Name	amylin, partial
NCBI Official Synonym Full Names	islet amyloid polypeptide
NCBI Official Symbol	IAPP
NCBI Official Synonym Symbols	DAP; IAP
NCBI Protein Information	islet amyloid polypeptide
UniProt Protein Name	Islet amyloid polypeptide
UniProt Gene Name	IAPP
UniProt Synonym Gene Names	DAP
UniProt Entry Name	IAPP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/amylin/180758>