

Canine Hyaluronidase ELISA Kit | HAase elisa kit

Canine hyaluronidase (HAase) ELISA Kit

CATALOG #
AAA180778

REACTIVITY
Canine


PRODUCT SYNONYMS

Hyaluronidase; Canine hyaluronidase (HAase) ELISA Kit; HAase elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	463
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HAase elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine HAase monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	2935328
NCBI GeneID	8372
Molecular Weight	15,865 Da

NCBI Official Full Name	hyaluronidase
NCBI Official Synonym Full Names	hyaluronoglucosaminidase 3
NCBI Official Symbol	HYAL3
NCBI Official Synonym Symbols	LUCA3; HYAL-3; LUCA-3
NCBI Protein Information	hyaluronidase-3
UniProt Protein Name	Hyaluronidase-3
UniProt Gene Name	HYAL3
UniProt Synonym Gene Names	LUCA3; Hyal-3; LuCa-3
UniProt Entry Name	HYAL3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/haase-canine-elisa-kits/hyaluronidase/180778>