

Canine Tissue Inhibitors Of Metalloproteinase 1 ELISA Kit | TIMP-1 elisa kit

Canine Tissue Inhibitors Of Metalloproteinase 1 (TIMP-1) ELISA Kit

CATALOG #
AAA180792

REACTIVITY
Canine


PRODUCT SYNONYMS

Tissue Inhibitors Of Metalloproteinase 1; Canine Tissue Inhibitors Of Metalloproteinase 1 (TIMP-1) ELISA Kit; TIMP-1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Sequence Length	40
Samples	Canine Serum, Plasma Or Cell Culture Supernatant And Organizations In The Natural And Recombinant Tc Concentration
Assay Type	Sandwich
Detection Range	20 mmol/L-0.312 mmol/L
Sensitivity	0.06 mmol/L.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TIMP-1 elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Canine TC monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

994731

NCBI GeneID	7076
Molecular Weight	23,171 Da
NCBI Official Full Name	tissue inhibitor of metalloproteinase 1, partial
NCBI Official Synonym Full Names	TIMP metallopeptidase inhibitor 1
NCBI Official Symbol	TIMP1
NCBI Official Synonym Symbols	EPA; EPO; HCI; CLGI; TIMP; TIMP-1
NCBI Protein Information	metalloproteinase inhibitor 1
UniProt Protein Name	Metalloproteinase inhibitor 1
UniProt Gene Name	TIMP1
UniProt Synonym Gene Names	CLGI; TIMP; EPA; Collagenase inhibitor; TIMP-1
UniProt Entry Name	TIMP1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/timp-1-canine-elisa-kits/tissue-inhibitors-of-metalloproteinase-1/180792>