

Canine Interleukin 6 ELISA Kit | IL-6R elisa kit

Canine Interleukin 6 receptor (IL-6R) ELISA Kit

CATALOG #
AAA180839

REACTIVITY
Canine


PRODUCT SYNONYMS

Interleukin 6; Canine Interleukin 6 receptor (IL-6R) ELISA Kit; IL-6R elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	212
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IL-6R elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine IL-6R monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	10834984
NCBI GeneID	3569
NCBI Accession #	NP_000591.1

NCBI GenBank Nucleotide #	NM_000600.4
Molecular Weight	23,718 Da
NCBI Official Full Name	interleukin-6 isoform 1
NCBI Official Synonym Full Names	interleukin 6
NCBI Official Symbol	IL6
NCBI Official Synonym Symbols	CDF; HGF; HSF; BSF2; IL-6; BSF-2; IFNB2; IFN-beta-2
NCBI Protein Information	interleukin-6
UniProt Protein Name	Interleukin-6
UniProt Gene Name	IL6
UniProt Synonym Gene Names	IFNB2; IL-6; BSF-2; CDF; IFN-beta-2
UniProt Entry Name	IL6_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/il-6r-canine-elisa-kits/interleukin-6/180839>