

Canine Angiotensinogen ELISA Kit | AGT elisa kit

Canine Angiotensinogen (AGT) ELISA Kit

CATALOG #
AAA180915

REACTIVITY
Canine


PRODUCT SYNONYMS

Angiotensinogen; Canine Angiotensinogen (AGT) ELISA Kit; AGT elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	485
Assay Type	Sandwich
Detection Range	100 ng/mL-1.56 ng/mL
Sensitivity	Up to 0.5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AGT elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine AGT monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	532198
NCBI GeneID	183

Molecular Weight	53,154 Da
NCBI Official Full Name	angiotensinogen
NCBI Official Synonym Full Names	angiotensinogen
NCBI Official Symbol	AGT
NCBI Official Synonym Symbols	ANHU; SERPINA8
NCBI Protein Information	angiotensinogen
UniProt Protein Name	Angiotensinogen
UniProt Gene Name	AGT
UniProt Synonym Gene Names	SERPINA8; Ang I; Ang II; Ang III; Ang IV
UniProt Entry Name	ANGT_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/agt-canine-elisa-kits/angiotensinogen/180915>