

Canine Thioredoxin ELISA Kit | Trx elisa kit

Canine Thioredoxin (Trx) ELISA Kit

CATALOG #
AAA180936

REACTIVITY
Canine


PRODUCT SYNONYMS

Thioredoxin; Canine Thioredoxin (Trx) ELISA Kit; Trx elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Sequence Length	105
Samples	Canine serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	2000 pg/mL-31.2 pg/mL
Sensitivity	Up to 12 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Trx elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine Trx monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

[16758644](#)

NCBI GeneID	116484
NCBI Accession #	NP_446252.1
NCBI GenBank Nucleotide #	NM_053800.3
Molecular Weight	11,673 Da
NCBI Official Full Name	thioredoxin
NCBI Official Synonym Full Names	thioredoxin 1
NCBI Official Symbol	Txn1
NCBI Official Synonym Symbols	Txn
NCBI Protein Information	thioredoxin
UniProt Protein Name	Thioredoxin
UniProt Gene Name	Txn
UniProt Synonym Gene Names	Txn1; Trx
UniProt Entry Name	THIO_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/trx-canine-elisa-kits/thioredoxin/180936>