

Canine Apolipoprotein B100 ELISA Kit | Apo-B100 elisa kit

Canine Apolipoprotein B100 (Apo-B100) ELISA Kit

CATALOG #
AAA180953

REACTIVITY
Canine


PRODUCT SYNONYMS

Apolipoprotein B100; Canine Apolipoprotein B100 (Apo-B100) ELISA Kit; Apo-B100 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	4563
Assay Type	Sandwich
Detection Range	1000 ng/mL-15.6 ng/mL
Sensitivity	Up to 5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Apo-B100 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine Apo-B100 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	105990532
NCBI GeneID	338

NCBI Accession #	NP_000375.2
NCBI GenBank Nucleotide #	NM_000384.2
Molecular Weight	515,605 Da
NCBI Official Full Name	apolipoprotein B-100
NCBI Official Synonym Full Names	apolipoprotein B
NCBI Official Symbol	APOB
NCBI Official Synonym Symbols	FLDB; LDLCQ4; apoB-48; apoB-100
NCBI Protein Information	apolipoprotein B-100
UniProt Protein Name	Apolipoprotein B-100
UniProt Gene Name	APOB
UniProt Synonym Gene Names	Apo B-100; Apo B-48
UniProt Entry Name	APOB_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/apo-b100-canine-elisa-kits/apolipoprotein-b100/180953>