

Canine anti-sperm antibody IgG (AsAb IgG) ELISA Kit | AsAb IgG elisa kit

Canine anti-sperm antibody IgG (AsAb IgG) ELISA Kit

CATALOG #
AAA181020

REACTIVITY
Canine


PRODUCT SYNONYMS

anti-sperm antibody IgG (AsAb IgG); Canine anti-sperm antibody IgG (AsAb IgG) ELISA Kit; AsAb IgG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Samples	Serum, plasma or cell culture supernatant and organizations in the natural and recombinant AsAb IgG concentration
Assay Type	Sandwich
Detection Range	20 ng/ml-0.312ng/ml
Sensitivity	0.06 ng/ml.
Intra-assay Precision	? 8%
Inter-assay Precision	? 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AsAb IgG elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antigen is Canine AsAb IgG antigen .Samples and antigen are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.