

Canine Galectin 3 ELISA Kit | GAL-3 elisa kit

Canine Galectin 3 (GAL-3) ELISA Kit

CATALOG #
AAA181125

REACTIVITY
Canine


PRODUCT SYNONYMS

Galectin 3; Canine Galectin 3 (GAL-3) ELISA Kit; GAL-3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	250
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GAL-3 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine GAL-3 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	115430223
NCBI GeneID	3958
NCBI Accession #	NP_002297.2

NCBI GenBank Nucleotide #	NM_002306.3
Molecular Weight	26,152 Da
NCBI Official Full Name	galectin-3 isoform 1
NCBI Official Synonym Full Names	lectin, galactoside-binding, soluble, 3
NCBI Official Symbol	LGALS3
NCBI Official Synonym Symbols	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2
NCBI Protein Information	galectin-3
UniProt Protein Name	Galectin-3
UniProt Gene Name	LGALS3
UniProt Synonym Gene Names	MAC2; Gal-3; CBP 35; GALBP
UniProt Entry Name	LEG3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/gal-3-canine-elisa-kits/galectin-3/181125>