

Bovine Complement 3b (C3b) ELISA Kit | C3b elisa kit

Bovine Complement 3b (C3b) ELISA Kit

CATALOG #
AAA181218

REACTIVITY
Bovine


PRODUCT SYNONYMS

Complement 3b (C3b); Bovine Complement 3b (C3b) ELISA Kit; C3b elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Specificity	No cross-reaction with other factors.
Sequence Length	1663
Samples	Serum, plasma or cell culture supernatant and organizations in the natural and recombinant C3b concentration
Assay Type	Sandwich
Detection Range	100 ng/ml-1.56 ng/ml
Sensitivity	0.5 ng/ml.
Intra-assay Precision	? 8%
Inter-assay Precision	? 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for C3b elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Bovine C3b monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	115298678
NCBI GeneID	718

NCBI Accession #	NP_000055.2
NCBI GenBank Nucleotide #	NM_000064.3
Molecular Weight	187,148 Da
NCBI Official Full Name	complement C3 preproprotein
NCBI Official Synonym Full Names	complement C3
NCBI Official Symbol	C3
NCBI Official Synonym Symbols	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p
NCBI Protein Information	complement C3
UniProt Protein Name	Complement C3
UniProt Gene Name	C3
UniProt Synonym Gene Names	CPAMD1; C3bc; ASP
UniProt Entry Name	CO3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/complement-3b-c3b/181218>