

Human Inducible T-Cell Co Stimulator (ICOS) ELISA Kit | ICOS elisa kit

Human Inducible T-Cell Co Stimulator (ICOS) ELISA Kit

CATALOG #
AAA183948

REACTIVITY
Human


PRODUCT SYNONYMS

Inducible T-Cell Co Stimulator (ICOS); Human Inducible T-Cell Co Stimulator (ICOS) ELISA Kit; CD278; AILIM; Inducible T-Cell Costimulator; Activation-Inducible Lymphocyte Immunomediatory Molecule; ICOS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Inducible T-Cell Co Stimulator (ICOS). No significant cross-reactivity or interference between Inducible T-Cell Co Stimulator (ICOS) and analogues was observed.
Sequence Length	199
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Urine, Cell Lysates and Other Biological Fluids
Detection Range	15.6-1,000pg/mL
Sensitivity	< 5.8pg/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Inducible T-Cell Co Stimulator (ICOS) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Inducible T-Cell Co Stimulator (ICOS) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ICOS elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Inducible T-Cell Co Stimulator (ICOS). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Inducible T-Cell Co Stimulator (ICOS). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Inducible T-Cell Co Stimulator (ICOS), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Inducible T-Cell Co Stimulator (ICOS) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	9968296
NCBI GeneID	29851
UniProt Accession #	Q9Y6W8; Q8N6W8
Molecular Weight	19,127 Da
NCBI Official Full Name	inducible T-cell co-stimulator
NCBI Official Synonym Full Names	inducible T cell costimulator
NCBI Official Symbol	ICOS
NCBI Official Synonym Symbols	AILIM; CD278; CVID1
NCBI Protein Information	inducible T-cell costimulator
UniProt Protein Name	Inducible T-cell costimulator
UniProt Gene Name	ICOS
UniProt Synonym Gene Names	AILIM

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.