

Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit | ABCG2 elisa kit

Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit

CATALOG #
AAA184213

REACTIVITY
Mouse


PRODUCT SYNONYMS

ATP Binding Cassette Transporter G2 (ABCG2); Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit; CDw338; CD338; ABC-G2; MRX; BCRP1; ABC15; ABCP; BCRP; BMDP; MXR; MXR1; Mitoxantrone resistance-associated protein; Breast Cancer Resistance Protein 1; Urate exporter; ABCG2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of ATP Binding Cassette Transporter G2 (ABCG2). No significant cross-reactivity or interference between ATP Binding Cassette Transporter G2 (ABCG2) and analogues was observed.
Sequence Length	611
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.055ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level ATP Binding Cassette Transporter G2 (ABCG2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level ATP Binding Cassette Transporter G2 (ABCG2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ABCG2 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to ATP Binding Cassette Transporter G2 (ABCG2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to ATP Binding Cassette Transporter G2 (ABCG2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain ATP Binding Cassette Transporter G2 (ABCG2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of ATP Binding Cassette Transporter G2 (ABCG2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	383792176
NCBI GeneID	9429
NCBI Accession #	NP_001244315.1
NCBI GenBank Nucleotide #	NM_001257386.2
UniProt Accession #	Q9UNQ0; O95374; Q4W5I3; Q53ZQ1; Q569L4; Q5YLG4; Q86V64; Q8IX16; Q96LD6; Q96TA8; A0A1W3; A8K1T5
Molecular Weight	67,453 Da
NCBI Official Full Name	ATP-binding cassette sub-family G member 2 isoform 2
NCBI Official Synonym Full Names	ATP binding cassette subfamily G member 2 (Junior blood group)
NCBI Official Symbol	ABCG2
NCBI Official Synonym Symbols	MRX; MXR; ABCP; BCRP; BMDP; MXR1; ABC15; BCRP1; CD338; GOUT1; MXR-1; CDw338; UAQTL1; EST157481
NCBI Protein Information	ATP-binding cassette sub-family G member 2
UniProt Protein Name	ATP-binding cassette sub-family G member 2
UniProt Gene Name	ABCG2
UniProt Synonym Gene Names	ABCP; BCRP; BCRP1; MXR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.