

Human Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) ELISA Kit | ENTPD1 elisa kit

Human Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) ELISA Kit

CATALOG # AAA184370	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1); Human Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) ELISA Kit; CD39; ATPDase; NTPDase-1; Ecto-ATP diphosphohydrolase 1; Ecto-apyrase; Lymphoid cell activation antigen; ENTPD1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1). No significant cross-reactivity or interference between Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) and analogues was observed.
Sequence Length	511
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.062ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ENTPD1 elisa kit

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Ectonucleoside Triphosphate Diphosphohydrolase 1 (ENTPD1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	12018242
NCBI GeneID	64519
NCBI Accession #	NP_072109.1
NCBI GenBank Nucleotide #	NM_022587.1
UniProt Accession #	P97687
Molecular Weight	57,408 Da
NCBI Official Full Name	ectonucleoside triphosphate diphosphohydrolase 1
NCBI Official Synonym Full Names	ectonucleoside triphosphate diphosphohydrolase 1
NCBI Official Symbol	Entpd1
NCBI Official Synonym Symbols	NTPDase-1
NCBI Protein Information	ectonucleoside triphosphate diphosphohydrolase 1
UniProt Protein Name	Ectonucleoside triphosphate diphosphohydrolase 1
UniProt Gene Name	Entpd1
UniProt Synonym Gene Names	Cd39; NTPDase 1; Ecto-ATPDase 1; Ecto-ATPase 1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.