

Human Fibroblast Activation Protein Alpha (FAPa) ELISA Kit | FAPa elisa kit

Human Fibroblast Activation Protein Alpha (FAPa) ELISA Kit

CATALOG #
AAA185748

REACTIVITY
Human


PRODUCT SYNONYMS

Fibroblast Activation Protein Alpha (FAPa); Human Fibroblast Activation Protein Alpha (FAPa) ELISA Kit; APCE; SIMP; DPPIV; Seprase; 170 kDa melanoma membrane-bound gelatinase; Post-proline cleaving enzyme; Serine integral membrane protease; Surface-expressed protease; FAPa elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Fibroblast Activation Protein Alpha (FAPa). No significant cross-reactivity or interference between Fibroblast Activation Protein Alpha (FAPa) and analogues was observed.
Sequence Length	239
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.781-50ng/mL
Sensitivity	< 0.31ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Fibroblast Activation Protein Alpha (FAPa) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Fibroblast Activation Protein Alpha (FAPa) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FAPa elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Fibroblast Activation Protein Alpha (FAPa). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Fibroblast Activation Protein Alpha (FAPa). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Fibroblast Activation Protein Alpha (FAPa), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Fibroblast Activation Protein Alpha (FAPa) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	119631759
NCBI GeneID	2191
UniProt Accession #	Q12884; O00199; Q53TP5; Q86Z29; Q99998; Q9UID4
Molecular Weight	26,954 Da
NCBI Official Full Name	fibroblast activation protein, alpha, isoform CRA_b
NCBI Official Synonym Full Names	fibroblast activation protein alpha
NCBI Official Symbol	FAP
NCBI Official Synonym Symbols	FAPA; SIMP; DPPIV; FAPalpha
NCBI Protein Information	prolyl endopeptidase FAP
UniProt Protein Name	Prolyl endopeptidase FAP
UniProt Gene Name	FAP
UniProt Synonym Gene Names	APCE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/fibroblast-activation-protein-alpha-fapa/185748>