

Human Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) ELISA Kit | CAMK2b elisa kit

Human Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) ELISA Kit

CATALOG #
AAA186227

REACTIVITY
Human


PRODUCT SYNONYMS

Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b); Human Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) ELISA Kit; CAM2; CAMK2; CAMKB; Proline Rich Calmodulin-Dependent Protein Kinase; CAMK2b elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b). No significant cross-reactivity or interference between Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) and analogues was observed.
Sequence Length	517
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.312-20ng/mL
Sensitivity	< 0.113ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CAMK2b elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Calcium/Calmodulin Dependent Protein Kinase II Beta (CAMK2b) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	26051218
NCBI GeneID	816
NCBI Accession #	NP_742081.1
NCBI GenBank Nucleotide #	NM_172084.2
UniProt Accession #	Q13554; O95437; O95438; O95599; Q9UGH7; A4D2K0; A4D2K1; A4D2K2; A4D2K3; A4D2K4; A4D2K5; A4D2K6
Molecular Weight	57,938 Da
NCBI Official Full Name	calcium/calmodulin-dependent protein kinase type II subunit beta isoform 8
NCBI Official Synonym Full Names	calcium/calmodulin dependent protein kinase II beta
NCBI Official Symbol	CAMK2B
NCBI Official Synonym Symbols	CAM2; CAMK2; CAMKB; MRD54
NCBI Protein Information	calcium/calmodulin-dependent protein kinase type II subunit beta
UniProt Protein Name	Calcium/calmodulin-dependent protein kinase type II subunit beta
UniProt Gene Name	CAMK2B
UniProt Synonym Gene Names	CAM2; CAMK2; CAMKB; CaM kinase II subunit beta; CaMK-II subunit beta

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.