

CLIP-associating protein 2 (CLASP2) Recombinant Protein | CLASP2 recombinant protein

Recombinant Human CLIP-associating protein 2 (CLASP2) , partial

CATALOG # AAA18663	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
------------------------------	---


PRODUCT SYNONYMS

CLIP-associating protein 2 (CLASP2); Recombinant Human CLIP-associating protein 2 (CLASP2); partial; CLASP2 recombinant protein

PRODUCT SPECIFICATIONS


Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Length	1,515
Species	Human
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

IMAGES




SDS-PAGE

Human CLIP-associating protein 2(CLASP2) ,partial His-tagged Yeast derived



SDS-PAGE

Human CLIP-associating protein 2(CLASP2) ,partial His-tagged Mammalian cell derived



SDS-PAGE

Human CLIP-associating protein 2 (CLASP2) ,partial His tag and fused with a special protein E.coli derived

NCBI AND UNIPROT INFORMATION

NCBI GI #	333440451
NCBI GeneID	23122

NCBI Accession #	NP_001193973.1
NCBI GenBank Nucleotide #	NM_001207044.1
UniProt Accession #	O75122 ; Q7L8F6 ; Q8N6R6 ; Q9BQT3 ; Q9BQT4 ; Q9H7A3 ; Q9NSZ2 ; B2RTR1 ; F5H604
Molecular Weight	165,924 Da
NCBI Official Full Name	CLIP-associating protein 2 isoform 2
NCBI Official Synonym Full Names	cytoplasmic linker associated protein 2
NCBI Official Symbol	CLASP2
NCBI Protein Information	CLIP-associating protein 2
UniProt Protein Name	CLIP-associating protein 2
UniProt Gene Name	CLASP2
UniProt Synonym Gene Names	KIAA0627; hOrbit2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/clasp2-recombinant-protein-e-coli-or-yeast-or-baculovirus-or-mammalian-cell/18663>