

Mouse Death Associated Protein Kinase 3 (DA PK3) ELISA Kit | DAPK3 elisa kit

Mouse Death Associated Protein Kinase 3 (DA PK3) ELISA Kit

CATALOG #
AAA186747

REACTIVITY
Mouse


PRODUCT SYNONYMS

Death Associated Protein Kinase 3 (DA PK3); Mouse Death Associated Protein Kinase 3 (DA PK3) ELISA Kit; ZIP; ZIPK; Dlk; DAP-like kinase; MYPT1 kinase; Zipper-interacting protein kinase; DAPK3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Death Associated Protein Kinase 3 (DA PK3). No significant cross-reactivity or interference between Death Associated Protein Kinase 3 (DA PK3) and analogues was observed.
Sequence Length	322
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.058ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Death Associated Protein Kinase 3 (DA PK3) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Death Associated Protein Kinase 3 (DA PK3) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DAPK3 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Death Associated Protein Kinase 3 (DAPK3). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Death Associated Protein Kinase 3 (DAPK3). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Death Associated Protein Kinase 3 (DAPK3), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Death Associated Protein Kinase 3 (DAPK3) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	4557511
NCBI GeneID	1613
NCBI Accession #	NP_001339.1
NCBI GenBank Nucleotide #	NM_001348.2
UniProt Accession #	O43293; Q05JY4; A0AVN4; B3KQE2
Molecular Weight	37,048 Da
NCBI Official Full Name	death-associated protein kinase 3
NCBI Official Synonym Full Names	death associated protein kinase 3
NCBI Official Symbol	DAPK3
NCBI Official Synonym Symbols	DLK; ZIP; ZIPK
NCBI Protein Information	death-associated protein kinase 3
UniProt Protein Name	Death-associated protein kinase 3
UniProt Gene Name	DAPK3
UniProt Synonym Gene Names	ZIPK; DAP kinase 3; Dlk; ZIP-kinase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/death-associated-protein-kinase-3-dapk3/186747>