

Human Calpain, Small Subunit 1 (CAPNS1) ELISA Kit | CAPNS1 elisa kit

Human Calpain, Small Subunit 1 (CAPNS1) ELISA Kit

CATALOG #
AAA186774

REACTIVITY
Human


PRODUCT SYNONYMS

Calpain; Small Subunit 1 (CAPNS1); Human Calpain; Small Subunit 1 (CAPNS1) ELISA Kit; CANP; 30K; Calpain 4; CANPS; CAPN4; CDPS; Calcium-activated neutral proteinase small subunit; Calcium-dependent protease small subunit; Calpain regulatory subunit; CAPNS1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Calpain, Small Subunit 1 (CAPNS1). No significant cross-reactivity or interference between Calpain, Small Subunit 1 (CAPNS1) and analogues was observed.
Sequence Length	268
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates, Cell Culture Supernates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.066ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Calpain, Small Subunit 1 (CAPNS1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Calpain, Small Subunit 1 (CAPNS1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CAPNS1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Calpain, Small Subunit 1 (CAPNS1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Calpain, Small Subunit 1 (CAPNS1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Calpain, Small Subunit 1 (CAPNS1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Calpain, Small Subunit 1 (CAPNS1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502565
NCBI GeneID	826
NCBI Accession #	NP_001740.1
NCBI GenBank Nucleotide #	NM_001749.3
UniProt Accession #	P04632; Q8WTX3; Q96EW0; A8K0P1
Molecular Weight	28,316 Da
NCBI Official Full Name	calpain small subunit 1 isoform 1
NCBI Official Synonym Full Names	calpain small subunit 1
NCBI Official Symbol	CAPNS1
NCBI Official Synonym Symbols	CANP; CDPS; CSS1; CANPS; CAPN4; CALPAIN4
NCBI Protein Information	calpain small subunit 1
UniProt Protein Name	Calpain small subunit 1
UniProt Gene Name	CAPNS1
UniProt Synonym Gene Names	CAPN4; CAPNS; CSS1; CANP small subunit; CDPS

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.